

//VMD 8 February 2021

//Multipoint OPD Counter Toggles

//Allows you to toggle between different counter types in the multipoint tool using shortcut keys in order to make distinctions between different microstructural features faster and easier when assessing osteon population densities

```
macro "Intact Osteon... [1]" {  
    run("Point Tool...", "type=Dot color=Yellow size=Medium label show counter=0");  
}
```

```
macro "Fragmentary Osteon... [2]" {  
    run("Point Tool...", "type=Dot color=Magenta size=Medium label show counter=1");  
}
```

```
macro "Howships Lacuna... [3]" {  
    run("Point Tool...", "type=Dot color=Cyan size=Medium label show counter=2");  
}
```

```
macro "Drifting Osteon... [4]" {  
    run("Point Tool...", "type=Dot color=Orange size=Medium label show counter=3");  
}
```

```
macro "Other... [5]" {  
    run("Point Tool...", "type=Dot color=Green size=Medium label show counter=4");  
}
```